





Aerospray® Gram Automated Slide Stainer/Cyocentrifuge (Model 7322)

Applications Training


Gram Stain Overview

- ◆ Developed ~125 years ago by Hans Christian Gram.
- ◆ Still widely used throughout microbiology.
- ◆ Stains bacteria, parasites, yeast (fungi), etc.
- ◆ Differentiates bacteria into two large groups, which is used as a first level classification.
 - ◆ Gram Positive (blue/purple)
 - ◆ Gram Negative (pink/red)
- ◆ Performed on body fluids or biopsies when infection is suspected.




Preparing a Gram Stain

- ◆ Cell concentration
 - ◆ Conventional centrifugation or
 - ◆ Cytocentrifugation (Cytopro Rotor)
- ◆ Fix Slides
 - ◆ Heat
 - ◆ Alcohol (methanol or ethanol)
- ◆ Apply Crystal Violet reagent
- ◆ Apply Mordant (Iodine), precipitates with CV in the cells
- ◆ Wash out precipitate with alcohol. CV precipitate washes out slower in Gram positive bacteria because of the thick peptidoglycan cell membrane.
- ◆ Counter stain Gram negative bacteria pink.



Preparing a Gram Stain


- ◆ Fixation (optional) – Methanol (PUMP E)
- ◆ Primary Stain – Crystal Violet (PUMP C)
- ◆ Wash – H₂O (PUMP D)
- ◆ Mordant – Iodine (PUMP B)
- ◆ Wash – H₂O (PUMP D)
- ◆ Decolorization/Counterstain (PUMP A)
- ◆ Wash – H₂O (PUMP D)
- ◆ Spin Dry



When to Perform a Gram Stain?



- ◆ Preliminary classification of causative agent to initiate a fast presumptive diagnosis and treatment
- ◆ Antibiotic / Antimycotic treatment monitoring
- ◆ Pre-screening of sputum samples
- ◆ Identification of bacteria in sweat (Johnson & Johnson)
- ◆ Screening of water lines for beverage suppliers (Coca Cola, Gatorade)




Contaminated sputum mixed flora (no predominance).





Advantages of Cytopro Rotor with Gram Stainer

- ◆ Concentrating samples with Cytopro is optional, but recommended when possible.
- ◆ Concentration of body fluids for superior sensitivity
- ◆ Uniform thin slide for reproducible staining result
- ◆ Fast (save 20 minutes/sample)



Gram Staining Reagents

◆ Reagent A: Decolorizer/Counterstain


- Decolorizer and Counterstain is added together. Options:
 - Safraine dye: mostly used in US and Asia.
 - ✓ Safranin Dye with Methanol and Isopropanol (SS-041A)
 - ✓ Safranin Dye with Isopropanol and Acetone (SS-041AA) – EBS preferred.
 - Fuchsin dye: mostly used in Europe.
 - ✓ Fuchsin Dye with Methanol and Isopropanol (SS-041AF)
 - ✓ Fuchsin Dye with Isopropanol and Acetone (SS-041AAF)

◆ Reagent B: Mordant- Iodine


- 6 month shelf life for ready to use stain (SS-041B)
- 8 month shelf life for concentrate (SS-141B)

◆ Reagent C: Primary stain- Crystal violet


E.Coli -safranin



E.Coli -fuchsin



S.aureus-fuchsin



Gram Staining Reagents




- ◆ There are four different decolorizer/counterstain reagents that can be used.

Decolorizer	Description
SS-041A– Safranin	<ul style="list-style-type: none">• Including 60% Isopropanol, 40% Methanol• More orange G-
SS-041AA – Acetone with Safranin	<ul style="list-style-type: none">• Including 75% Isopropanol, 25% Acetone• Gives stronger gram positive result.• EBS recommended
SS-041AF - Fuchsin	<ul style="list-style-type: none">• Including 60% Isopropanol, 40% Methanol• Used mainly in Europe• More purple G-
SS-041AAF – Acetone with Fuchsin	<ul style="list-style-type: none">• Including 75% Isopropanol, 25% Acetone• Gives stronger gram positive results.



7322 Prepared & Concentrate Reagents



Gram Staining Reagents



1 Place the reagent bottles front to back in the following order:

- (A) Decolorizer with Counterstain
- (B) Iodine
- (C) Crystal Violet
- (D) Deionized Water
- (E) Methanol or Ethanol

Programming Flexibility



- ◆ 9 decolorizer settings.
- ◆ Crystal violet (CV) and Iodine (I2) adjustments
 - ◆ High settings of CV and I2 will make Gram positive bacteria stronger (more positive)

STAIN PROGRAM SETTINGS

Program Name: Example Program


Decolorizer: 9

Fixation: Off Normal High

Crystal Violet: Low Medium High

Iodine: Low Medium High

Save



Programming Flexibility









- ◆ Most customers will find a setting that works well for all of their samples.

Decolorizer Setting	Thickness	Suggested Specimen
1	Very thin Smears	Old culture(s)
2	Thin Smears	CSF, urine, peritoneal, vaginal, wounds, bronchial washes, etc.
3	Thin and medium smears	Vaginal, sputum, bronchial washed, CSF, urine, wounds, etc.
4		
5	Medium and Thick Smears	Sputum, bronchial washes, blood cultures, etc.
6		
7	Thick smears	Sputum, bronchial washes, blood cultures, etc.
8		
9	Very Thick Smears	Stools, blood cultures, etc.



7322 Decolorizer/ Counterstain Comparison



Quality Control Management



ELITechGroup Gram Control slides SS-250 (50 slides)

- ◆ For quality control management
- ◆ Gram positive: S.aureus – Gram negative: E.coli



Advantages of Gram Stainer

- ◆ No cross contamination
- ◆ Fast, easy, CLEAN
- ◆ Consistent staining
- ◆ Economical, more slides with less stain.
 - Much less stain waste generated.
- ◆ 12 slide or 30 slide carousel
- ◆ Stainer/Cytocentrifuge combination
- ◆ Flexibility
 - Nine different decolorizer settings
 - Increased Crystal Violet and/or Iodine spray
 - Optional fixation
- ◆ Minimal maintenance



7322 Disadvantages



- ◆ **Over/under decolorization**
 - ◆ Over decolorization: Gram positive organisms stain pink
 - ◆ Under decolorization: Gram negative organisms stain purple
- ◆ **Nozzle maintenance: Iodine (B line) needs to be checked or cleaned once a week.**
- ◆ **Bowl fouling with Crystal Violet and Iodine precipitate**
- ◆ **Short Expiration date of Iodine (six months)**
- ◆ **Monthly B-line flush**




Overcoming the 7322 Disadvantages




◆ Over/under decolorization

- ◆ Decolorizer setting of 3 or 4 are comparable to hand staining with 30 or 40 seconds of stationary decolorizer treatment.
- ◆ It is well known that gram positive cultures, such as *Bacillus* sp. And *Streptococcus* sp., can easily be over-decolorized, especially when cultures are old or stressed.
- ◆ The 7322 is more sensitive to these cultures than hand staining. (Remind customers to use fresh cultures (log phase) when these specimens are used as controls)
- ◆ Samples at the upper and lower edges of the slide may over-decolorize because specimens do not get best stain coverage.



Competition: Hand Stains

- ◆ Inconsistent results
- ◆ More reagent consumption
- ◆ Time consuming (separate decolorizer/counterstain step + slides not dry)
- ◆ Contact with hazardous stain
- ◆ Not clean, messy sinks.
- ◆ No traceability and not user-friendly



Competition: PREVI Color Gram

- ◆ EBS is manufacturer
- ◆ EBS has much more knowledge & service available.
- ◆ More expensive (device+reagents) as barcode reader, 5L reagent level and 10L waste level detection is not optional
- ◆ Samples can only be entered via barcode, not manually
- ◆ 6 nozzle starter kit, EBS 2 nozzle starter kit
- ◆ BioMérieux cannot offer full stainer range, only Gram Stainer and TB with different technology



Competition: Dip Stainers

◆ Dip stainers in general:

- Cross-contamination (batch for 60 samples)
- More reagent consumption
 - Dispose reagents or filter them if cell washoff is expected.
- Reagents can be exposed to environment.
- Fumes from reagents.
- Wet slides.
- Limited traceability.

◆ RAL

- \$0,50-\$1,00/slide (EBS 0,03\$/slide)
- Slower TAT (15' 40 slides – ELITechGroup 9' for 60 slides)
- Short shelf life of reagents once opened.
- Reagent waste.



Market Development

- ◆ Support is where we can differentiate ourselves.
- ◆ Evaluate collaborations with local suppliers (BD, Siemens) to offer a more complete solution in tenders.
- ◆ Focus on labs preferring manual sample entry.
- ◆ Always promote Cytopro in parallel as a niche solution!



Market Development

◆ Do you also collaborate with another partner?


- Latin American
 - Aerospray Gram + Phoenix (BD)
- Germany
 - Aerospray Gram + WASP® (Copan)
 - Total of 102 stainer in Germany, some placed along with WASP.



Questions to ask during sales visit



- ◆ How do you Gram stain your microbiology samples?
- ◆ What type of reagents does your lab use? (explain safranine/fuchsin, acetone or methanol)
- ◆ How many samples does your lab stain/day? (explain reagent calculation tool website – economical benefits)
- ◆ Do you sometimes encounter problems with uniformity of slides or cross-contamination? (explain Aerospray fully automated, standardized system)
- ◆ How do you treat your waste? (explain limited reagent consumption)
- ◆ Does your internal QC management or accreditation process require traceability for patients and reagents.
- ◆ How do you concentrate your body fluids? (explain Cytopro option)
- ◆ Would you be interested in a dual purpose cost-efficient staining/cytocentrifugation system?
- ◆ How is re-imburement arranged (per test or global budget)?



Performing a Successful Demonstration



◆ Prospection:


- Understand the customer's fixation protocol and preferred staining reagents.
- Show EBS stained control slides to preliminary select preferred decolorizer setting.
- Present Cytopro rotor for standardized sample preparation

◆ Demo preparation:

- Ask customer to prepare fresh samples (max 24h) . Record 'abnormal' characteristics before staining.
- Include SS-250 Control slides.

◆ Actual demo:

- Start with Decolorizer setting 4 or setting customer preferred on stained control slide.
- Staining control slides (SS-250), can help determine preferred Decolorizer setting.
- For new accounts, set up stainer with safranine + acetone (SS-041AA).
- Demo on Cytopro rotor.




Performing a Successful Demonstration




◆ Case study 1: Question

- ◆ The **Jyntendo Hospital in Japan** evaluated the Gram stainer. They used acetone with safranin and obtained the following stain result at decolorizer setting 7.



- ◆ How does the slide look?
- ◆ What can cause the slide to look like this?
- ◆ How can we modify Aerospray Gram settings to meet the customer need?



Performing a Successful Demonstration



◆ Case study 1: Answer

- Over-decolorization
- Solutions:
 - First exclude error on test operation, but all ok
 - Was all methanol removed from lines during set-up?
 - Reagents correctly prepared?
 - Nozzles clean?
 - Reagent grade acetone used?
 - Fresh sample, lag phase?
 - Decolorizer setting should be decreased to setting 3
 - In addition, Crystal Violet/Iodine can be increased to « High »

Note: G+ (Bacillus, Strep) are more sensitive to over-decolorization




Performing a Successful Demonstration



◆ Case study 2: Question

- ◆ The Policia Lab in Colombia evaluated the Gram stainer. They used safranin with acetone. They observed over-decolorization and debris on slide.



- ◆ What can cause debris to be seen on slide and how can we prevent it?




Performing a Successful Demonstration



◆ Case study 2: Answer

- Causes: expired iodine used
- Other possible causes:
 - B line (iodine) not properly cleaned (monthly flush). Iodine debris can be present on the slide.
 - D nozzles (water) not properly cleaned.
 - ◆ D line bottle should be cleaned with bleach (sodium hypochloride) one per month. If not cleaned the line can be contaminated causing debris on the slide.
 - Dirty slides used. Debris stained on the slide.
 - Fresh sample?




Performing a Successful Demonstration




◆ Case study 3: Question

- ◆ The ARUP Lab in Salt Lake City Utah evaluated the Gram stainer. They observed under-decolorization on a control slide. (Purple stain on a lot of G- organisms.



- ◆ What can cause under-decolorization? How is it corrected.



Performing a Successful Demonstration



◆ Case study 3: Answer

- ◆ Causes: Decolorizer setting is too low. Increase decolorizer setting number.
- ◆ Other possible causes:
 - ◆ Check A nozzle, make sure it is spraying properly.

