



ELITechGroup

EMPOWERING IVD



Biomedical Systems

MACRODUCT® ADVANCED

THE ONLY CYSTIC FIBROSIS SWEAT TEST FULLY COMPLIANT TO GUIDELINES



Macroduct® Advanced

SWEAT COLLECTION SYSTEM

- 1 Touchscreen user interface with built-in step-by-step graphical instructions eliminating user-error.
- 2 Revised electrodes, collector, and straps for greater convenience when working with small limbs.
- 3 Collector is clearly marked for easy determination of minimum sweat volume and sufficient sweat rate.
- 4 The perfect companion for CF diagnostic of babies and adults.



Chlorocheck®

SWEAT CHLORIDE ANALYZER

MACRODUCTADVANCED.COM