



## AEROSPRAY® Control Slides

Providing quality control to your laboratory


Innovation with Integrity



# Valuable tip for your quality management

Standardized control slides are very useful for establishing new methods in the laboratory. Furthermore, control slides are very useful for routine quality control (recommended for every run), protocol validation and assay optimization.

**Can be used to verify quality control management techniques.**  
**Support staff training and use of new reagents.**  
**No need to stock cultures.**



Control slides  
(50 slides/package)

## TB Quality control slides

Catalog #	Description
SS-249	TB stain control slides, pack of 50

Positive control: **Inactivated *Mycobacterium tuberculosis***  
Negative control: ***Escherichia coli***

## Gram Quality control slides

Catalog #	Description
SS-250	Gram stain control slides, pack of 50

Positive control: **Gram positive *Staphylococcus aureus***  
Negative control: ***Escherichia coli***

Ten of the corresponding control slides are included with the purchase of Aerospray® Gram or TB stainer.




**ELITechGroup Inc.**  
Logan, UT · USA

1.800.453.2725

© 2025 ELITechGroup, Inc. All rights reserved.  
Aerospray® is a registered trademark of ELITechGroup Inc.



**Online information**  
[Aerospraystaining.com](https://www.aerospraystaining.com)



BMS-60-0080-01A

